



Cellulose Acetate Membrane Electrophoresis Tank



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Horizontal Electrophoresis

Darsun Scientific

Cellulose Acetate Membrane Electrophoresis Tank

DS-CAME-7



Product Features

DS-CAME-7 is upgraded from single-bridge Double-bridge design.
Can keep a constant temperature with a cooling device.
The removable rod can be adjusted to meet different distance. Removable electrodes enable maintained and replacement easy and quick.
Auto-switch-off when the lid is moved

MODEL	DS-CAME-7
Gel size (WxL)	85x200 (mm)
No. of samples	10~364
Total buffer volume	~50 (ml)
Dimension (Lx W x H)	300x150x61 (mm)
Weight	0.5 (kg)

DS-CAME-7 Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
85x200mm Gel tray	1
Baffle	2
Combined comb frame	1
1.0mm thickness 10-wells comb	12
1.8mm thickness 18-wells comb	12