



Conductivity Meter



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Conductivity Meter

DDS-TZSA



Conductivity Meter

Model		DDS-TZSA
Parameters		Conductivity/TDS/°C
Range	Cond	1.0 μ S/cm – 100 mS/cm
	TDS	0 – 1000 mg/L
	Temp	(0.0 ~ 60.0) °C
Accuracy	Cond	±0.5%FS
	TDS	±1.0%FS
	Temp	±0.3 °C ±1bit
Temp. Compensation		(0.0 ~ 40.0) °C, auto/manual
Power		220V / 110V
Dimension / N.W.		300 × 200 × 90 mm / 1.5 kg
Shipping size / G.W.		330 × 260 × 280 mm / 2 kg
Standard Accessories		Electrode Holder
		Cond. electrode (K=1)
		Power adapter
		Manual
		Temperature probe