



Conductivity Meter



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Conductivity Meter

DSS-ELA



Conductivity Meter

Model		DSS-ELA
Parameters		Conductivity
Range	Cond	(0 ~ 2 × 10 ⁵) μS/cm
	TDS	—
	Temp	—
Accuracy	Cond	±1.0%FS ±1bit
	TDS	—
	Temp	—
Temp. Compensation		(5 ~ 35) °C, manual
Power		220V / 110V
Dimension / N.W.		300 × 200 × 90 mm / 1.5 kg
Shipping size / G.W.		330 × 260 × 280 mm / 2 kg
Standard Accessories		Electrode Holder
		Cond. electrode (K=1)
		Power adapter
		Manual
		—