




Cellulose Acetate Membrane Electrophoresis Tank



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Cellulose Acetate Membrane Electrophoresis Tank

DS-CAME-7A



Product Overview

Can keep a constant temperature with a cooling device.
The removable rod can be adjusted to meet different distance. Removable electrodes enable maintained and replacement easy and quick.
Auto-switch-off when the lid is moved.

Model	DS-SP7A
Gel size (WxL)	66x250 (mm)X2
	82X250 (mm)X2
Total buffer volume	~560 (ml)
Dimension (Lx W x H)	318x175x 65 (mm)
Weight	1.0 (kg)

Partsname	Quantity
Lid with cables	1
Body tank with electrodes	1
Adjusting rod	4
Cellulose acetate Membrane	10
Cooling tube	2