



DS-Mini Transfer (Classical)



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

DS-Mini Transfer (Classical)

DS-ZYZ10T



Product Features

Simultaneously transfer 1-2 precast or hand-cast gels within an hourThe blue ice box is completely contained within this KIT and absorbs heat generated during quick transfersThe tank of the Tranfer Cell can also be used as the tank of our Vertical cell Perfect match with DS Vertical cell KIT (DS-DYC10T)

Product Overview

The DS Transfer Cell (DS-ZYZ10T) is a part of the Protein Electrophoresis system. It works with the vertical cell (DS- DYC10V) together to provide high-quality blotting of the hand-cast and precast gels. This model accommodates two gel holder cassettes for the transfer of the gels run in the vertical cell (DS-DYC10V).

Model	DS-ZYZ10T
Size	75x100mm
Time	about60min
Qty of the gels	1-2
Cooling method	Ice Box
Buffer Volume	450ml