



HM

Heating Magnetic Stirrer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Heating Magnetic Stirrer

DSHS-5



FEATURES



Strong magnetic field



Stainless steel ceramic coated countertops

Model	DSHS-5
Control Mode	Brushless motor
Max Stirring Volume	2L(H ₂ O)
RPM	200-1500rpm
Timing Range	/
Control Mode	Knob
Working Plate (mm)	Φ135
Plate Temp.	260℃
Allowable Ambient Temp	5~40℃ 80%RH
Power	200W
Product Size (mm)	225*155*90