

Horizontal Electrophoresis Tank



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific




Horizontal Electrophoresis

Darsun Scientific

Horizontal Electrophoresis Tank

DS-SHET-G



Product Features

Abundant buffer not only guarantees cooling effect s, but also keeps the pH value stable during the whole experiment process.
Removable electrodes enable maintained and replacement easy and quick.
Adding sample more quick and suitable for the multi- channel pipettes.
With a special gel casting, preventing the gel leakage, operation is simple and convenient.
It's more convenient to open the lid with the easy-open buttons.
It prevents buffer volatilizing and electricity leakage with transparent lid.
Auto-switch-off when the lid is moved.

MODEL	DS-SHET-G
Gel size (WxL)	245x250 (mm)
No. of samples	19~364
Total buffer volume	~200 (ml)
Dimension (Lx W x H)	487x300x93 (mm)
Weight	3.3 (kg)

DS-SHET-G Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
245x250mm Gel tray	1
Gel casting stand	1
1.0mm thickness 19-wells comb	2
1.8mm thickness 19-wells comb	2
1.0mm thickness 26-wells comb	2
1.8mm thickness 26 wells comb	2
1.0mm thickness 44 wells comb	2
1.8mm thickness 44 wells comb	2
1.0mm thickness 52- wells comb	2
1.8mm thickness 50-wells comb	2
Level adjusting button	4