



LC

Large Capacity Magnetic Stirrer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Large Capacity Magnetic Stirrer

DSMS-5L/DSMS-10L
DSMS-15L/DSMS-30L
DSMS-50L



FEATURES



Timing function



Knob operation



Forward and reverse



Large capacity

Model	DSMS-5L	DSMS-10L	DSMS-15L	DSMS-30L	DSMS-50L
Control Mode	Display & Knob				
Display Method	LCD display				
Power	10W	20W	20W	30W	50W
Working Plate (mm)	180*180	220*220	220*220	260*260	330*330
Max Stirring Volume	5L	10L	15L	30L	50L
Speed Range	50-1500rpm	50-1500rpm	50-1500rpm	50-1300rpm	50-1100rpm
Speed Regulation Accuracy	5rpm				
Timing Range	1min-23h59min				
Stir Bar Size (mm)	40	50	60	80	100
Product Size (mm)	280*180*65	315*220*65	315*220*65	350*260*85	435*330*85