



OS

Overhead Stirrer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Overhead Stirrer

DSOS-5/DSOS-10/DSOS-15S
/DSOS-30S



FEATURES



Heat more uniformly



304 stainless steel plate,
Corrosion resistant



LCD display



Durable brushless motor

Model	DSOS-5	DSOS-10	DSOS-15S	DSOS-30S
Motor	Brushless motor			
Max Amount Of Stirring	5L	10L	10L	20L
Range Of Rotate Speed	200-1500rpm	100 - 1000rpm	100 - 2000rpm	50 - 1200rpm
Timing	/	/	1min-23h59min	1min-23h59min
Max Torque	10N.cm	20N.cm	150N.cm	30N.cm
Max Viscosity	5000mPas	10000mPas	8000mPas	16000mPas
Drill Chuck Diameter	1-10mm			
Power	20W	25W	30W	36W
Voltage	AC100-240V			