



PB

Precision balance




India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Precision balance

DS-203EL-PLUS



Product Features

- . Internal auto calibration/strain gauge sensor
- . LCD(white back light with black font)
- . Tare function/counting/unit conversion
- . percentage/sensitivity/speed set
- . density/peak holding
- . print function (baud rate setting continues communication or press communication)
- . calibration time setting
- . overload alarm/fault alarm/level indicator
- . option:-USB/printer/density kit

Model	DS-203EL-PLUS
Capacity	220 g
Readability	0.001 g
Min weighing	0.004 g
Repeatability	±0.002 g
Linearity	±0.003 g
Stable Time	≤3 s
Operation temp	5–25 °C
Pan size	Φ80 mm (Standard) Φ90 mm (option)
Dimension	305×210×305 mm
Draft shield size	175×165×200 mm (Inner diameter)
LCD Size	94×24 mm
Cal. weight	Internal calibration
G.W	6.0 kg
N.W	4.5 kg
Packing size	425×310×410 mm
Interface	RS232 (Standard)
Sensor	Load cell