



Semi-dry Transfer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Horizontal Electrophoresis

Darsun Scientific

Semi-dry Transfer

DS-TB4



Product Features

Platinum-coated titanium alloy plate for the anode and stainless-steel plate for the cathode, which ensure reliable transfer effect.
Flat-plate anode with spring device, which can automatically adjust the pressure according to different thickness of gel, convenient to operation..
Simple hanging buckle design, which ensures that the stacked transfer "sandwiches" are completely connected with the electrodes without moving.
Multiple gels placed side by side or stacked can be transferred at one time. Fast speed, transfers can be completed within 60mins.
Minimal buffer required to soak the filter papers. Automatic power-off when the lid is opened.

MODEL	DS-TB-4
Trans size (WxL)	260x210(mm)
Dimension (LxWxH)	350x240x83 (mm)
Weight	5.2 (kg)