



Transilluminator



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Transilluminator

DS-02 -UV



Product Overview

No need for darkroom, can be used under different weather.
Cut-gel windows for cutting the gel.
UV filter glass, quartz ultraviolet lamp, long life service.

Model	DS-02 UV
Transmission size (Wx L)	200x150 (mm)
Transmission UV wavelength	302 (nm)
Reflectance UV wavelength	254x365 (nm)
Transmission UV lamp power	302nm (8W)
Reflectance UV lamp power	254nm (11W) 365nm (11W)
Dimension (LxWxH)	435x295x490 (mm)
Weight	16.9 (kg)

Parts name	Quantity
Mainframe	1
Camera rack	1
Power connecting cables	3
5A insurance tube	2
Power cord	1