



HS

Digital Display Hotplate Stirrer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Digital Display Hotplate Stirrer

DSHS-5C/DSHS-5P



FEATURES



Strong magnetic field



Stainless steel ceramic coated countertops

Model	DSHS-5C	DSHS-5P
Control Mode	Brushless motor	Brushless motor
Max Stirring Volume	2L(H ₂ O)	2L(H ₂ O)
RPM	200-1500rpm	50-1500rpm
Timing Range	/	1min-12h
Control Mode	Digital display & Knob	Digital display & Button
Working Plate (mm)	Φ135	Φ135
Plate Temp.	280°C	310°C
Liquid/Plate Temp Accuracy	±1°C/±3°C	±1°C/±3°C
Power	300W	400W
Product Size (mm)	225*155*90	225*155*90