

Horizontal Electrophoresis Tank



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632


DarsunScientific

Horizontal Electrophoresis

Darsun Scientific

Horizontal Electrophoresis Tank

DS-SHET-C



Product Features

Abundant buffer not only guarantees cooling effects, but also keeps the pH value stable during the whole experiment process. Removable electrodes enable maintained and replacement easy and quick. With the fluorescent ruler on the tray, can be observed directly into the ultraviolet. With a special gel casting, preventing the gel leakage, operation is simple and convenient. It's more convenient to open the lid with the easy-open buttons. It prevents buffer volatilizing and electricity leakage with transparent lid. Auto-switch-off when the lid is moved.

MODEL	DS-SHET-C
Gel size (WxL)	130x130 (mm)
	78x100 (mm)
	48x75 (mm)
No. of samples	3~104
Total buffer volume	~600 (ml)
Dimension (Lx W x H)	305x152x80 (mm)
Weight	1.7 (kg)

DS-SHET-C Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
130x130mm Gel tray	1
78x100mm Gel tray	1
48x75mm Geltray	2
Gel casting stand	1
1.0mm thickness 13-wells comb	2
1.5mm thickness 13-wells comb	2
1.0mm thickness 26-wells comb	4
1.0mm thickness 26 wells comb	2
1.0mm thickness 10 wells comb	2
1.0mm thickness 15 wells comb	2
1.5mm thickness 3- wells comb	2
1.0mm thickness 5-wells comb	2
Level adjusting button	4