

## Horizontal Electrophoresis Tank

---



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific




# Horizontal Electrophoresis

Darsun Scientific

## Horizontal Electrophoresis Tank

DS-SHET-F



### Product Features

Abundant buffer not only guarantees cooling effect s, but also keeps the pH value stable during the whole experiment process.  
Removable electrodes enable maintained and replacement easy and quick.  
Adding sample more quick and suitable for the multi- channel pipettes.  
With a special gel casting, preventing the gel leakage, operation is simple and convenient.  
It's more convenient to open the lid with the easy-open buttons.  
It prevents buffer volatilizing and electricity leakage with transparent lid.  
Auto-switch-off when the lid is moved.

MODEL	DS-SHET-F
Gel size (WxL)	200x200 (mm)
	200x150 (mm)
	200x100 (mm )
No. of samples	17~104
Total buffer volume	~800 (ml)
Dimension (Lx W x H)	397x230x93 (mm)
Weight	3.3 (kg)

### DS-SHET-F Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
200x200mm Gel tray	1
200x150mm Gel tray	1
200x100mm Geltray	1
Gel casting stand	1
1.0mm thickness 17-wells comb	2
1.5mm thickness 17-wells comb	2
1.0mm thickness 22-wells comb	2
1.0mm thickness 22 wells comb	2
1.0mm thickness 36 wells comb	2
1.0mm thickness 30 wells comb	2
1.5mm thickness 44- wells comb	2
1.0mm thickness 44-wells comb	2
Level adjusting button	4