

Horizontal Electrophoresis Tank



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific




Horizontal Electrophoresis

Darsun Scientific

Horizontal Electrophoresis Tank

DS-SHE-3



Product Features

Fast electrophoresis running speed, the whole process can be finished within 15mins. 96 sample detection in one time. Be quicker to add samples and suitable for multi-channel pipettes. It's more convenient to open the lid with the easy-open buttons. Auto-switch-off when the lid is moved.

MODEL	DS-SHE-3
Gel size (WxL)	85x200 (mm)
No. of samples	10-364
Total buffer volume	~50 (ml)
Dimension (Lx W x H)	300x150x61 (mm)
Weight	0.5 (kg)

DS-SHE-3 Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
85x200mm Gel tray	1
Baffle	2
Combined comb frame	1
1.0mm thickness 10-wells comb	12
1.8mm thickness 18-wells comb	12