



LS

Low Speed Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



FEATURES



Temperature protection
& over shock protection



Double electronic
lock



Program memory



Multi rotors
optional

Darsun Scientific

Low Speed Centrifuge



Model	DSLC500-8	DSLC500-12	DSLC500-18	DSLC500-24	DSLC500-6
Rotor Capacity	15ml*8 Centrifuge tube	15ml*12 Centrifuge tube	15/50ml*6 Centrifuge tube	2/5/10ml*18	2/5/10ml*24
	2/5/10ml*8 Blood tube	2/5/10ml*12 Blood tube	2/5/10ml*6 Blood tube	Blood tube	Blood tube
Speed	500 - 4000rpm				
Max RCF	1967xg	2075xg	2146xg	2182xg	2164xg
Semidiameter Of Rotor	11cm	11.6 cm	12cm	12.2 cm	12.1 cm
Timing Range	30s - 59min50s				
Protection Function	Lid opening protection, cavity over temperature protection, vibration protection				
Power	60W	60W	80W	80W	80W
Power Supply	AC 110V - 220V, 50/60Hz				