

## Vertical Electrophoresis Tank



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific




# Horizontal Electrophoresis

Darsun Scientific

# Vertical Electrophoresis Tank

DS-VE5/ DS-VE5L



### Product Features

The installation of glass plates can be completed within 15s, easy and quick. With the gel casting modular multi plates can be made . It's more convenient to open the lid with the easy-open buttons. Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process. Auto-switch-off when the lid is moved.

MODEL	DS-VE5	DS-VE5L
Glass plate size (WxL)	216X110(mm)	216x220 (mm)
Gel size (WxL)	186x95(mm)	186x205 (mm)
Gel thickness	1.0 (mm)	1.0(mm)
Number of gel	1~2 pc (s)	1~2 pc (s)
No of sample	25~104	25~104
Total buffer volume	~2200(ml)	~8000 (ml)
Dimension (LxWxH)	300X190X129 (mm)	320x240x230 (mm)
weight	3.2 (kg)	5.6 (kg)

Lid with cables	1
Body tank	1
Gel casting core with electrodes	1
Gel casting modular	2
Notched glass plate	4
Spacer plate with 1.0mm glass spacer	4
1.0mm thickness 25-will comb	2
1.0mm thickness 40-will comb	2
1.0mm thickness 52-will comb	2
Dropper	1