



# Analyzer DS-181



## India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam  
Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

## Australia Headquarters

21b Lawler Drive, Oran Park,  
NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific

Darsun Scientific

# Analyzer

DS-181



## Product Features

Product Type: Bench-Top Fully automatic chemiluminescence immunoassay Analyzer

Product weight: Net weight: 62 kg

Processor: RAM 2GB, main frequency 4 × 1 GHz or above

Security Interface: USB 2.0 & Ethernet RJ45 interface

1. Total 90 test parameters including infectious markers (HIV, HCV, HBV & Syphilis), Dengue range (NS1 Ag, IgM & IgG), cardiac markers, fertility markers, etc.
2. Solid Phase Carrier: Magnetic Particles.
3. Processing mode: Batch & continuous random access.
4. Result Throughput is 120 tests/hour (max.), first result within 15 min.
5. 30 sample position & 10 reagent positions to process multiple test samples at a time.
6. Patented Single steel S/R Probe (needle) for sample & reagent addition contributing in high reliability & more accuracy in test results.
7. Single assay cup for reagent & sample processing leading to lower cost in consumables.
8. Large data storage capacity for storing around 5, 00,000 patient data.
9. HIS & LIS compatible with application support.

**DARSUN**  
SCIENTIFIC

**Darsun**Scientific