



Constant Temperature Circulator



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Constant Temperature Circulator

DS-1050



Product Features

DS-1050 Constant Temperature Circulator used for control keeping experiment's temperature at stable values during electrophoresis run s, DS-1050 also can be connected with the tank or other lab equipment to serve as a constant temperature controller.

Model	DS-1050
Temperature control range	-5~35 (°C)
Temperature accuracy	+ 0.1 (°C)
Circulating fluidrate (L/ min)	5 (L/MIN)
Tank volume (ml)	2.2 (L)
Power	300 (W)
Dimension (LxWxH)	196x406x482 (mm)
Weight	20.0 (kg)

Parts name	Quantity
Mainframe	1
Power cord	1