



Gel Vacuum Drier



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Gel Vacuum Drier

DS-GV-A



Product Features

DS-GV-A Gel Vacuum Drier is designed for all kinds of electrophoresis gel drying. In order to prevent gel cracking , gel drier is used for processing it during daily experiments.

Model	DS-GV-A
Drying size (WxL)	Drying size (WxL)
Drying control	30°C-80°C
Timing range	20~120min
Volume (LxWxH)	556x380x60 (mm)
Weight	15.5 (kg)

Parts name	Quantity
Mainframe	1
Filter paper	1
5A insurance tube	2
Powercord	1