



LED TRANSILLUMINATOR



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

LED
TRANSILLUMINATOR

DS-01 UV-II



Model	DS-01 UV-II
Transmission size (WxI)	200x50 (mm)
Transmission UV wavelength	254x365 (nm)
Transmission UV lamp power	8 (W)
Dimension (LxWxH)	270x78x110 (mm)
Weight	0.55