



WB

Water Bath



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park, NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific



Darsun Scientific

Water Bath

DS-1-WB100/DS-2-WB100
DS-4-WB100/DS-6-WB100
DS-8F-WB100



FEATURES



One-piece stainless steel tank



High corrosion resistance



High precise temperature control



Digital display
thermostatic

Model	DS-1-WB100	DS-2-WB100	DS-4-WB100	DS-6-WB100	DS-8F-WB100
Temperature Range	Room temp-100°C				
Power	300W	600W	800W	1200W	1500W
Accuracy	±0.1°C				
Working Chamber Size (mm)	160*160*130	300*150*140	300*300*140	450*300*130	620*300*120