



# HARD TOP CHEST FREEZER

---



## India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

## Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific

# Darsun Scientific HARD TOP CHEST FREEZER

DSC-400DM5/DSC-588DM5/DSC-588DM5



## Product Feature

These 2-in-1 hard top chest freezers from Darsun conveniently double up as chillers at the turn of a knob and are available in a wide variety of sizes ranging from 145 to 780 litres to suit every conceivable need. You can also choose between single and double lid variants.

Model	DSC-400DM5	DSC-500DM5	DSC-700DM5
Star Rating	5	5	5
Gross Volume (Litres)	383	524	734
Dimensions wxdxh (mm)	1235x745x905	1645x745x905	1860x855x905
Temperature Range	≤-20°C	≤-20°C	≤-20°C
Power (W):	300	380	300
Refrigerant:	R290	R290	R290
Lock	Yes	Yes	Yes
Defrost	Manual	Manual	Manual
No. of Lids / Doors	2	2	2
No. of Baskets	1	1	1
Wheels / Castors	Yes	Yes	Yes
Certification	CE/ISO		