



Manual colony counter



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam
Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

Australia Headquarters

21b Lawler Drive, Oran Park,
NSW,
2570, Australia
info@darsunscientific.com
+61-283135514, 430583632

DarsunScientific

Darsun Scientific

Manual colony counter

DS-A11



Product Overview

1. Displays average values
2. Bright LED lighting
3. Three-digit display for easy readings Adjustable pressure
4. Beeping noise prevents errors

Product Features

Manual colony counter DS-A11 accommodates a wide range of colony counts from 0 to 999 for accurate measurements. This colony counter fits petri dishes from 50 to 150 mm, offering versatility. It features multi-mode functionality for simultaneous counting from two petri dishes. Equipped with a pressure-sensitive touch-tip pen marker. Ensures thorough examination and accurate counting of colonies through cumulative marking. Provides nine times magnifications for easier counting.

Model	DS-A11
Counting range	0 to 999
Compatible petri dish	50 to 150 mm
Magnification	3 to 9 times
Digital display	3 digits
Rack	2
Relative Moisture	0.8
Ambient temperature	5 to 50°C
Protective class	IP21
Power consumption	20 W
Power	AC 100 to 240V/ 50 to 60 Hz
Dimension	360x300x180mm
Weight	4.0 kg