



**SD**

## Spray Dryer-5L



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific



Darsun Scientific

Spray Dryer-5L

DSF-5L



- The drying speed is rapid. After centrifugal spraying, the surface area of the liquid is greatly increased, and 95-98% of the water can be evaporated in a moment in the high temperature air flow, and the drying time needs more than ten seconds.
- Parallel flow spray drying is adopted. During the drying process, the temperature of the droplets is not high, and the product quality is good.
- Due to the instantaneous completion of the drying process, the particles of the product can basically maintain a spherical shape similar to liquid droplets, and the product has good dispersibility, fluidity, and solubility.
- Simplified production process and convenient operation control.

Model	DSF-5L
Air inlet temperature control	30~300°C
Air outlet temperature contra	30~140°C
Maximum feed amount	Max. 5000ml/h (adjustable)
Nozzle speed	36000r/h (adjustable)
Feeding mode	Peristaltic pump adjustment
Minimum feed amount	50ml
Average drying time	1.0~1.5S
Draught fan	3kW
Electric Heater Power	9kW
Machine power	13KW /380V
Overall dimension	1600X1200X2000( mm )