

## Largest capacity refrigerated Centrifuge

---




### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific



# Darsun Scientific

## Largest capacity refrigerated Centrifuge

DSL-8M



### Product feature

1. Widely used in the field of blood bank, pharmaceutical factory and laboratory.
2. Brushless frequency motor in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, 10 kinds of program, 10 kinds of acceleration and deceleration for choice.
5. There are different rotors for choice.
6. Electric lid lock, compact design, super speed and imbalance protection.
7. The centrifuge body is made of high-quality steel, safe and reliable.

Model	DSL-8M
Parameter	Specification
Max Speed	8000rpm
Timer	1~9h59min
Speed Accuracy	±20rpm
Noise	≤65DB
Power supply	AC380V/220V 50HZ/60HZ
Max Volume	8x2000ml (16x500ml blood bag)
Dimension (L×W×H)	940x890x1000mm
Max RCF	14551xg
Net Weight	600kg
Acc/Dec	10 kinds
Temperature Range	-20°C~40°C
Temperature Accuracy	±1°C

Matched Rotors				
Order No	Rotor	Max Speed (r/min)	Volume (ml)	Max RCF (×g)
DS8M-1	Swing Rotor	4000	8x2000ml (16x500ml blood bags)	5716
DS8M-2	Swing Rotor	4000	6x2400ml	5330
DS8M-3	Angle Rotor	8000	6x1000ml	14551
DS8M-4	Swing Rotor	4200	2x6x1000ml	5680