



## VM

# Vortex Mixer



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### Australia Headquarters

21b Lawler Drive, Oran Park, NSW,  
2570, Australia  
[info@darsunscientific.com](mailto:info@darsunscientific.com)  
+61-283135514, 430583632

**Darsun**Scientific



# Darsun Scientific

# Vortex Mixer

DSM25/DSM25D



- Simple and reliable structure, small dimension, low power consumption and low noise;
- Two models optional: fixed and high/low speed;
- Used for various mixing application with optional adaptor;
- Accept a variety of test sample container: centrifuge tube/PCR plate/cell culture plate/deep well plate;
- Stable operation, small speed error and long service life.

Model	DSM25	DSM25D
Power	30W	30W
Load Adapter	Flat vibrating head	Flat vibrating head / Microplate clamp
Speed Range	2800 rpm	L: 1300 rpm / H: 2800 rpm
Vibration Amplitude	4 mm	
Mix Method	Orbital	
Working Mode	Continuous	
Dimensions (L×W×H)	160 × 160 × 115 mm	
Weight	3.75 kg	
Power Supply	110V/220V; 50/60Hz	