



OS

3D Oscillator Shaker



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

3D Oscillator Shaker

DSSK-3D



FEATURES



Timing function



Brushless motor with strong power



3D mixing



Suitable for a variety of containers

Model	DSSK-3D
Type	3D Oscillating
Display	LCD
Tit Angle	10°
Max Load	0.8kg
Speed	10 - 80rpm
Timing	1min - 8h
Allowed Ambient Humidity Temp.	5~40°C 80%RH
Voltage	AC 100-240V, 50/60Hz
Power	10W
Dimension (mm)	184*184*130