



Semi-Auto Biochemistry Analyzer,



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Semi-Auto Biochemistry Analyzer,

DS-100C



Product Features

1. 7" inch true color touch screen(800*480 pixels)
2. 5 USB ports & 1 ethernet port
3. Flash: 256 MB;Over 300 programs, 200000 results can be stored.
4. Resolution: 0.0001 Abs.
5. Photometric range: 0.0000-4.0000 Abs.


Product Overview

Semi-Auto Biochemistry Analyzer is suitable for the analysis of serum, plasma, urine and other body fluid samples, and can complete the test of biochemical, special protein, drug monitoring and blood coagulation items



Multifunctional Incubator With Accurate Constant Temperature

Multifunctional incubator specially designed for chemistry and coagulation tests, different holes can be used to put sample tubes, chemistry cuvettes, coagulation cuvettes, accurate temperature control, more efficient sample incubation



- 5 USB port & Ethernet port
- Mouse, keyboard and printer can be connected. freely easy transmission of data



Windows embedded compact 7 smart operation system

More convenient to edit test items freely, clear presentation of real-time curve, calibration, quality control and sample results.



7" inch true color touch screen
(800*480) pixels