



HA

Automatic 5-Diff Hematology Analyzer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific



Darsun Scientific

Automatic 5-Diff Hematology Analyzer

DSK-8520



Double Chambers for WBC and RBC/PLT
29 parameter + 4 histograms + 1 Scatter Diagram
Up to 90 Samples/hour

- Semiconductor laser scattering combined with cytochemical staining technology to achieve accurate five-type white blood cell classification
- Advanced dual-sheath flow device for precise counting and classification of cells
- Excellent abnormal cell alarm function, customizable alarm rules
- Floating cluster classification technology improves the ability to identify the boundary of the scatter plot and makes the white blood cell classification more accurate

Model	DSK-8520
Display	12.1 inches colorful touch-screen display
Storage	Up to 200000 results with graphs
Data Transportation	USBX4, RS232,LAN HDMI, Duplex LIS
Language	Multi language can be customized
Print	Multiple print formats and can be customized
Dimension	55X40X53 cm
Weight	31 kg