



AC

Automatic Chemistry Analyzer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States


DarsunScientific



Darsun Scientific

Automatic Chemistry Analyzer

DS-6090



- Discrete, random stat function
- ISE module (for option , K Na, Cl, Ca, PH)
- Up to 400 tests/hour ,
- 93 sample positions
- 40 R1 & R2 reagent position (with cooling system)
- 90 sample cuvettes with auto washed system
- 2 independent mixing stirrers
- Read spectrophotometry, 12 wavelengths
- Probe anti-collision prevention
- Easy to operate and ,maintain windows 10 compliant
- Dimension:- 95X71X118 cm
- Weight:- 136 kg