



HA

Automatic veterinary 3-diff Hematology Analyzer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific



Darsun Scientific

Automatic veterinary 3-diff Hematology Analyzer DSK-VET3



Model	DSK-VET3
21 parameters	WBC LYM#,MID#, GR# LYM% MID% ,GR% ,RBC,HGB,HCT,MCV,MCH RDW-CV, RDW-SD, PLT,MPV, PDW, PCT, P-LCC,P-LCR
3 Histograms	WBC,RBC,PLT
Species	12 preset species including Canine ,Equine , Bovine, Murine , Swine , Lipoid 4 optical species
Sample volume	Whole blood 10 µl
Storage capacity	Auto-save , more than 200.000 analysis results available
Throughout	60 tests/hour
Fate search mode	By sample ID, patient ID, date, multi-keywords. Loading/extracting data, offline reading
Display	10.4 inch colorful touch screen
Operation mode	L-J ,X-R,X-B, etc.
Printout	Thermal print (defaulted) available for external printers optional printing templates and format
Data communication	Wifi, Bluetooth
Ports	LAN,X,RS232X1,HDMIX1,USB port x4
Language	English
Optional Accessories	Handheld barcode scanner
Dimension & weights	400x350x440mm 17 kg
Power supply	110-240V , 50/60Hz