



DB

## Dry Bath Shaking Incubator



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States


**Darsun**Scientific



Darsun Scientific

# Dry Bath Shaking Incubator

DS-DB100C/DS-DB100



LCD display



Over temperature protection



Program memory, 9 temperature points



Two modules, freely combined

| Model                 | DS-DBS100C  | DS-DBS100    |
|-----------------------|---|--------------|
| Temperature Range     | 0°C ~ 100°C   | 15°C ~ 100°C |
| Temperature Stability | ±0.5°C  | ±0.5°C       |
| Display Accuracy      | 0.1°C   | 0.1°C        |
| Heating Rate          | 5.5°C/min   | 5.5°C/min    |
| Cooling Rate          | 2.5°C/min(100°C-Room temperature) 0.5°C/min(Below room temperature) | /            |
| Speed / Timing        | 200-1500rpm / 1s-23h59min   |              |
| Mixing Diameter       | 3mm   |              |
| Program Memory        | 9 grades  |              |
| Power                 | 200W  |              |
| Voltage               | 100-240V 50/60HZ  |              |
| Size                  | 200*235*120mm   |              |

| module options | Numbers and aperture of test tube   |
|----------------|---|
| DBs1           | 0.2ml*96 96-hole PCR microplate   |
| DBs2           | 1.5ml *24 centrifuge tube   |
| DBs3           | 2ml *24 centrifuge tube   |
| DBs4           | 15ml *8 centrifuge tube   |
| DBs5           | 50ml *4 centrifuge tube   |
| oem of modules | Modules can be customized according to customer needs   |
| DBS-T          | Low module support for thermal cover<br>(thermal cover and transparent house cannot be used at the same time) |