



## GEL DOCUMENTATION SYSTEM (CLASSICAL)

---



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

# Darsun Scientific

## GEL DOCUMENTATION SYSTEM (CLASSICAL)

### DSH-510 /DSH-520



#### Product Features

- UP to 20 Megapixels
- Auto Focus
- Real White LED transmission
- Numeric Focus and Aperture

#### SPECIFICATION

Model		DSH-510	DSH-520
Digital Camera	Resolution	5Megapixels, 2592*1944	20Megapixels, 4800*4440
	Dynamic Range	4.8OD	
	Pixel density	16bit (65536 Grey Scales)	
Lens	Motorized Lens	F1.2 Motorized Lens	
	Resolution	20 Megapixels	
	Control	Numeric Focus controlled via Software	
Light Source	UV transmission	302nm as standard, Image Area : 21x26cm (Optional 254nm,365nm)	
	White LED transmission	Real white LED transmission, Image Area:19x26cm	
	Optional Blue LED transmission	Optional Blue LED transmission, Image Area:19x26cm	
Filter system	Filter	Standard filter for gel doc application	
Software	Capture	Capture software	
	Analysis	Analysis software	