



# Gel Vacuum Drier



## India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

## USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

Darsun Scientific

# Gel Vacuum Drier

DS-GV-A



## Product Features

DS-GV-A Gel Vacuum Drier is designed for all kinds of electrophoresis gel drying. In order to prevent gel cracking, gel drier is used for processing it during daily experiments.

Model	DS-GV-A
Drying size (WxL)	Drying size (WxL)
Drying control	30°C-80°C
Timing range	20-120min
Volume (LxWxH)	556x380x60 (mm)
Weight	15.5 (kg)

Partsname	Quantity
Mainframe	1
Filter paper	1
5A insurance tube	2
Powercord	1

