



**MO**

## Microplate oscillator



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

## FEATURES



Timing function



Brushless motor  
with strong power

# Darsun Scientific Microplate oscillator

DSVM-500-1P  
DSVM-500-2P



Model	DSVM-500-1P	DSVM-500-2P
Motor Type	Brushless motor	
Horizontal Gyration Amplitude	4mm	
Speed	100-1200rpm	
RPM Increment	10rpm	
Timing	1min - 8h	
Handling range	ELISA Plate / Tissue Culture Plate / 96-well plate	
Allowed Ambient Humidity Temp	5 ~ 40°C 80%RH	
Voltage	AC 100 - 240V, 50/60Hz	
Power	20W	
Dimension (mm)	155*150*90	180*150*90