



Premium Vortex Mixer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Premium Vortex Mixer



Vortex Mixer provides comprehensive mixing capabilities for enzymatic and RIA assays, general test tube mixing, biochemical assays, viral dilution preparation, precipitated assays, cell suspension vortexing and tissue sample mixing.

Specification

- Continuous and Touch modes
- Speed Accuracy: $\pm 5\%$ rpm
- Convenient Speed adjustments by knob control
- Stable for long-term operations
- Multiple accessories for different vessels
- Digital Speed display
- Speed range: 200 to 3000rpm

Model	Premium Vortex Mixer
Speed	Up to 200–3000 rpm
Warranty	2 years
Load Bearing Capacity	0.5 kg
Shipping Weight	14.99 kg
Display	Digital
Dimensions (L x W x H)	290x270x240mm 11.4x10.7x9.5
Motor Type	DC motor
Hertz	50/60 Hz
Material	ABS housing
Includes	Rubber platform
Tube Capacity	≤ 19 tubes
Operating Modes	Touch, Continuous

