



# RM

## Roller mixer



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

## FEATURES



**Timing function**



**Anti - lock design, more convenient**



**Adjustable speed**



**Brushless motor with strong power**

# Darsun Scientific Roller mixer

DSRM-6H/DSRM-8H



Model	DSRM-6H	DSRM-8H
Motor Type	Brushless motor	
Display	LCD	
Type	Horizontal rolling	
Speed	20-80rpm	
Number Of Rollers	6	8
Roller Length (mm)	280	
Roller Center Distance(mm)	34	
Specs	φ 6-34mm	
Timing	1min-23h59min	
Allowed Ambient Humidity Temp.	5~40°C 80%RH	
Voltage	AC 100-200V,50?60Hz	
Power	25W	
Dimension (mm)	440*255*105	430*320*100