

Semi-dry Transfer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters


30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Semi-dry Transfer

DS-TB3



Product Overview

- Quick transferring it can be completed within 1 hour.
- Less heat, no need for extra cooling device.
- Less buffer is needed than the tank transfer.
- Auto-switch-off when the lid is moved.

MODEL	DS-TB3
Trans size (WxL)	140x140 (mm)
Dimension (LxWxH)	185x175x40 (mm)
Weight	1.4 (kg)