



Transilluminator



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Transilluminator

DS-02 -UV



Product Overview

No need for darkroom, can be used under different weather.

Cut-gel windows for cutting the gel.

UV filter glass, quartz ultraviolet lamp, long life service.

| Model | DS-02 UV |
|----------------------------|-------------------------|
| Transmission size (Wx L) | 200x150 (mm) |
| Transmission UV wavelength | 302 (nm) |
| Reflectance UV wavelength | 254x365 (nm) |
| Transmission UV lamp power | 302nm (8W) |
| Reflectance UV lamp power | 254nm (11W) 365nm (11W) |
| Dimension (LxWxH) | 435x295x490 (mm) |
| Weight | 16.9 (kg) |

| Parts name | Quantity |
|-------------------------|----------|
| Mainframe | 1 |
| Camera rack | 1 |
| Power connecting cables | 3 |
| 5A insurance tube | 2 |
| Power cord | 1 |