



Tests Fully Auto Chemistry Biochemistry Analyzer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

DS-200



1. Working Mode : End point, Fixed, Kinetic, Fixed type Kinetic
2. Filter : 7+ 3(Optional)
3. Flowcell : 30 ul quartz
4. Control: 3 Controls per test
5. Display: 640x480 apixel
6. Printer: Inbuilt
7. Memory: 100,000

Model NO.	DS-200
Group	All
Analysis Method	End-Point, Fixed-Time, Rate (Kinetic)
Reagent Positions	28 Reagent Positions
Reagent Cooling	Refrigerated Tray with Independent Switch
Mixer Probe	Independent Stirring
Spectrophotometry	Post-Spectral Spectrophotometry
Absorbance	0~3.0 ABS
Lis System	Bi-Direction, Support HI7 Protocol
Trademark	DARSUN
Production Capacity	100000 unit/Year
Certification	CE
Throughput	200 tests/Hour
Sample Positions	37 Sample Positions
Probe Washing	Automatic Washing Interior and Exterior
Cuvettes	48 Reusable Cuvettes, Optical Length 6mm
Light Source	6V/10W Halogen Lamp
Wavelength	340, 405, 450, 480, 505, 546, 570, 600, 630, 700, 750, 800 nm
Software	Windows 7/8/10, 32 or 64 Bit
Interface	LAN Port Access
Specification	Gross weight: 55 kg