



HS

# Digital Display Hotplate Stirrer



#### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

#### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

Darsun Scientific

# Digital Display Hotplate Stirrer

DSHS-8P



## FEATURES



Strong magnetic field



Stainless steel ceramic coated countertops

Model	DSHS-8P
Control Mode	Brushless motor
Max Stirring Volume	5L(H <sub>2</sub> O)
RPM	50-1500rpm
Timing Range	1min - 12h
Control Mode	Digital display & Knob & Button
Working Plate (mm)	Square Plate160*160mm
Plate Temp.	310°C
Allowable Ambient Temp	5~40°C 80%RH
Power	500W
Product Size (mm)	285*180*90