



HS High Speed Centrifuge



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

FEATURES



Temperature protection function



Enclosed high-speed type rotor



Program memory



Stable rotation

Darsun Scientific

High Speed Centrifuge

DSMC-13P DSMC-15P



Model	DSMC-13P	DSMC-15P
Control Mode	Digital display & Button	
Max Relative Centrifugal Force	11337xg	14104xg
Speed Range	1000-13000rpm±5%	1000-14500rpm±5%
Rotor Capacity	0.2/0.5/1.5/2ml*12 Centrifuge tube 0.5/1.0ml*12 Micro blood tube	
Max Speed Acceleration / Reduction	≤12 s/ ≤15s	
Noise Under Max Rotation Speed	≤47db (A)	
Input Power	50W	55W
Shell Protection Level	AC100-240V / DC24V	
Allow Surrounding Temperature	5 - 40°C / 80%RH	

DarsunScientific