



LS

# Low Speed Centrifuge



#### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

#### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

Darsun Scientific

# Low Speed Centrifuge

DSLCL-400



## FEATURES



**Electronic door lock**



**Program memory**



**Emergency switches  
for easier use**



**LCD display**

Model	DSLCL-400
Speed Range	500-4500rpm
Max Relative Centrifuge Force	2180xg
Rotor Capacity	15ml*4 Centrifuge tube 2/5/10ml*8 Blood tube
Increment	100rpm
Timing Range	30s ~ 59min50s
Noise	≤ 60dB(A)
Voltage	AC220V
Power	40W
Dimensions (mm)	265*260*190