



Continuous Sealing Machine



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific



Darsun Scientific

Continuous Sealing Machine

DS100-A



Model	DS100-A	
Technical Parameters	1.	Sealing Speed : 10±0.5m/min
	2.	Sealing Margin : 0-35mm adjustable
	3.	Sealing Line Width : 12mm
	4.	Sealing Strength Meet the requirements of YY/T 0698.5-2023
	5.	AC Power : 110V/220V±10% 50/60Hz
	6.	Maximum Power:500W
	7.	LED two-color high-light digital tube display, touch button
	8.	Automatic control, continuous sealing, temperature deviation less than 1%, working temperature 60-220°C arbitrary setting
	9.	Rapid heating design, room temperature rises to 180 °C less than 40 seconds, high efficiency and energy saving, to meet the needs of fast-paced work;
	10.	Advanced flat heating element, can dry burn, high temperature resistance, long service life, high thermal efficiency
	11.	The floating constant pressure pressing system is designed to meet the sealing requirements of paper-plastic bags, paper-plastic three-dimensional bags and paper-paper bags.
	12.	Positive sequence sealing counter, it can count the number of sealing packages within 0 -9999
	13.	Automatic fault alarm indication
	14.	Auxiliary cooling design, equipped with a microcomputer controlled cooling system, reduce the waiting time of temperature conversion from high temperature sealing to low temperature sealing;
	15.	High safety. When the deviation between the sealing temperature and the set working temperature exceeds 4 °C, the equipment will automatically stop working, effectively ensuring the sealing quality and safe operation of the equipment;
	16.	Housing Material Carbon steel (Default) / stainless steel / aluminium alloy
	17.	Indentation Style Flat wheel sealing (Default) / Single indentation sealing / Pattern wheel sealing
	18.	Optional Accessories Extended roller-able workbench, medical sealing machine workstation
	19.	Plug: Optional for European / American
Size (LxWxH)	490mmx240mmx120mm	
Weight N.W	11kg	