



Economical Freeze Dryer 12L-80



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific


Darsun Scientific

Economical Freeze Dryer

DS12L-80A DS12L-80B DS12L-80C
DS12L-80E DS12L-80D



Model	DS12L-80A DS12L-80B DS12L-80C DS12L-80D DS12L-80E
condenser temperature	-85°C
Ice condenser capacity	12L/24h
Dimension (mm)	518x860x707(without chamber)
Vacuum pressure (pa)	<1Pa (standard <5Pa)
vacuum pump	2XZ-2(2L/s)
*optional vacuum pump	DRV-10(10³/h), low noisy and better performance
Refrige system	2 compressors Cascade
Defrost	Heating coil
Rated power	2100w
Net weight (Appro)	130kg(without vacuum pump,vacuum pump 20kg)



More Optional



DRV vacuum pump

Ulvaco Vacuum pump

Oil free Vacuum pump



4 Trays + Acrylic chamber

3 stopper Assembly Trays and Acrylic Chamber



4 Trays + 8 Manifold valve Chamber

stopper Assembly Tray + Manifold valve Chamber

Ampule Bottle / Test tube Adaptor



India Office

16 S/F B/S, T-Block Extn.
 Jain Colony Part-2 Uttam Nagar,
 New Delhi-110059
 info@darsunscientific.in
 (+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
 Sheridan, Wyoming 82801
 United States