



Economical Freeze Dryer 4L-60



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Economical Freeze Dryer

DS4L-60A DS4L-60B DS4L-60C
DS4L-60E DS4L-60D



There is another appearance Pioneer 3L, with name Pioneer -3L -60-2, same specification only different appearance.



Model					
	DS4L-60A	DS4L-60B	DS4L-60C	DS4L-60D	DS4L-60E
condenser temperature	-60°C				
Ice condenser capacity	4L/24h				
Dimension (mm)	392X557X683 (without chamber)				
Vacuum pressure (pa)	<1Pa (standard <5Pa)				
vacuum pump	2XZ-2(2L/s)				
*optional vacuum pump	DRV-10(10 ³ /h)				
Refrige system	Auto Cascade				
Defrost	Heating coil				
Rated power	1350w				
Net weight (Appro)	55kg(without vacuum pump,vacuum pump 20kg)				

