



## Economical Freeze Dryer 4L-80

---



### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific

# Darsun Scientific

## Economical Freeze Dryer

DS4L-80A DS4L-80B DS4L-80C  
DS4L-80E DS4L-80D



Model					
	DS4L-80A	DS4L-80B	DS4L-80C	DS4L-80D	DS4L-80E
condenser temperature	-85°C				
Ice condenser capacity	4L/24h				
Dimension (mm)	518X860X707(without chamber)				
Vacuum pressure (pa)	<1Pa (standard <5Pa))				
vacuum pump	2XZ-2(2L/s)				
*optional vacuum pump	DRV-10(10³/h),low noisy and better performance				
Refrige system	2 compressors Cascade				
Defrost	Heating coil				
Rated power	1350w				
Net weight (Appro)	75kg(without vacuum pump , vacuum pump 20kg)				

