



Economical Freeze Dryer 6L-60



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

UsA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific

Economical Freeze Dryer

DS6L-60A DS6L-60B DS6L-60C
DS6L-60E DS6L-60D



Model					
	DS6L-60A	DS6L-60B	DS6L-60C	DS6L-60D	DS6L-60E
condenser temperature	-60°C(-80°C/-100°C check below)				
Ice condenser capacity	6L/24h				
Dimension (mm)	518X860X707(without chamber)				
Vacuum pressure (pa)	<1Pa (standard<5Pa))				
vacuum pump	2XZ-2(2L/s)				
*optional vacuum pump	DRV-10(10³/h),low noisy and better performance				
Refrige system	Auto Cascade				
Defrost	Heating coil				
Rated power	1200w				
Net weight (Appro)	(without vacuum pump,vacuum pump 20kg)				

